

#### Lower Dose Csl Scintillator | Fastest AED | 14" X 17" Cassette-Sized DR

The AirDR™ is equipped with a directly deposited Cesium iodide (Csl) scintillator which provides unparalleled image quality. The cassette-size AirDR ™minimizes patient dose and improves throughput, making it an innovative Fail-Safe Detector.

The AirDR™ is our most advanced, low-dose digital radiography solution when paired with the robust and feature-rich image processing XC™ software.

XC™ provides superior quality images every time with an option to further manipulate the images.

# iCRco | Image Capture Review

### IMAGE

#### AIRDR | 100 MICRON DR

AirDR provides the highest efficiency cesium-iodide based scintillator, 16 bit available gray scale, and market leading 100 micron active maxtrix that guarantees image resolution up to 5 lp/mm.

### CAPTURE

#### XC | ACQUISITION

XC touchscreen scquisition with ICE-4 Enhancement Processing provides all-new features including, "Image Display State" to ensure balanced presentation of both soft tissue, overlapping bone structures, and automatic analysis of image characteristics to optimize processing.

### REVIEW

#### **CLARITY PACS**

Our fully web-enabled and integrated PACS solution help transition your practice into a safe, secure, and filmless environment. Clarity PACS™ supports all your current and future imaging needs.

### AirDR Specifications



The AirDR System has been specially designed and optimized to advance the imaging equipment you're using right now. Utilizing its unique form-factor and embedded Automatic Exposure Detection (AED), the AirDR System is compatible with any X-ray system designed to work with ISO 4090 - compliant, 35 x 43 cm cassettes.

√ No need to modify your generator, or bucky

√ No need to replace your grids

√ No need to discard your wall stand, or table

√ Water and Dust Resistant IP-42 Rating √ Carbon Fiber Construction Durability

√ Same Day Installation

	Panel	Amorphous Silicon active TFT/diode array, Carbon-fiber construction
	Scintillator	Direct Deposit; Cesium Iodide
	Pixel Matrix	3556 x 4320
	Pixel Pitch	100 μm
	lmage Data	16 bit
lmag	ge Transfer Time	Wired: 500 ms; Wireless: 3000 ms
	Active Area	True 35.5 cm x 43.2 cm
Exte	rnal Dimensions	ISO 4090 cassette size 14" x 17" (38.4 cm (w) x 46 cm (l) x 1.5 cm (h))
	Weight	8.4 lbs (3.8 kg)
W	eight Capacity	Maximum load weight: 150 kg
St	tatus Display	LED display (Wifi/Battery/Sensor)
Wi	reless Data I/F	802.11n Wifi standard
M	/ired Data I/F	GigE via optional power & communication tether
	X-ray I/F	Automatic Exposure Detection (AED)
Limi	ting Resolution	5 lp/mm
	Typical MTF	70% (1 lp/mm), 40% (2 lp/mm), 15% (4 lp/mm) for RQAS
•	Typical DQE	75% (0 lp/mm), 60% (1 lp/mm), 40% (3 lp/mm) for RQAS
	nvironment	10 - 35 °C operating, 30 - 70 % RH operating (non-condensing)
	Battery	Rechargeable battery, 53.3 Wh
Ва	ttery Charger	External two bay charger 100-240 V AC 50/60 Hz
Interfa	ce and Power Unit	Optional AirDR IPU with external power supply 100-240 V AC, GigE, and X-ray I/F
	Standards	IEC 60601-1, IEC 60601-2, IEC 60601-1-6, FCC 47CFR PT 15, FCC OET 65C, ETSI EN 301 893, EN 62311, ISO 10993-5, ISO 10993-10, CE
Binne	ed Mode (option)	Up to 8 fps for 2 x 2 binned, 200 $\mu m$ pitch for a pixel matrix of 1778 x 2160"
Image	Calibration (option)	On-board offset, gain and defective pixel corrections
Fast	Preview (option)	4 x 4 binned quick preview image
© 2020 iCRco All rights reserved "True Flat Scan Path" and "YC' are registered trademarks of iCRco   RR 102816AUS   Design & Specifications		

© 2020 iCRco. All rights reserved. "True Flat Scan Path" and "XC' are registered trademarks of iCRco | BR 102816AUS | Design & Specifications are subject to change without notice.

## **XC<sup>TM</sup> Acquisition Software Features**

XC - Intuitive touchscreen acquisition

ICE- 4 Processing - Automated image characteristics analysis for maximum image enhancement

Image Display State: Automated enhancement of image display at the point of acquisition

Smart search, sort, and filter options

Integration with front office management systems like RIS and EMR

Full set of annotation and measurement tools

User-preferred settings and privileges

### Outer Dimensions

